

catoolRT 21.9 Release Notes

Changes since 21.6

- Support XCP v1.4 with DAQ_PACKED_MODE
- XCP use appropriate time units depending on expected rate
- Various XCP fixes including sending values as NaN if parameter no longer exists
- Add "Display Delay" and "Cycles Per Analysis" information fields
- Rename "Upload" to "Upload to Support"
- Better handle checking for NI hardware when it has recently been plugged in
- Support digital triggering for interpolated clocks. This enables interpolation of low-resolution encoders
- Correct anti-aliasing for frame rate font
- Fixes to support of digital input channels
- DAQ channel configuration window is now re-sizeable
- Show correct PFI/port channels for trigger input